

FAST ANALYSIS OF Ne, O₂, N₂, CH₄ AND CO

Technique : Micro-GC

Column : Tamis 5Å


Carrier gas : Argon

Column temperature : 130°C

Injection time : fixed

Column pressure : 22 psi

μTCD sensitivity : Standard



The Micro-GC is a powerful solution that provides fast, accurate and easy-to-use solution to analyze your gas sample.

A wide range of capillary columns and MEMS injectors are available to optimize channel combinations inside the Micro-GC to obtain the desired separation. The μTCD is the high-performance universal detector used in such systems. A parallel analytical configuration of the channels or analytical modules will allow the sampling of a small amount of gas sample and simultaneous injection in all channels to develop up to 4 chromatograms and one cumulative report..

SRA Instruments is the ideal Partner with specific competences to develop with you turn key solutions based on optimized Micro-GC, accessories and dedicated software utility or report calculations for specific plant context.

Compound	RT (sec)	Concentration (%)
Ne	15.73	0.043%
O ₂	18.15	4.941%
N ₂	20.13	5.334%
CH ₄	25.29	3.807%
CO	26.65	4.760%

Some examples of applications are:
 gas coke analysis,
 biogas, syngas,
 refinery and fuel gas,
 natural gas, fuel cells.



More information : www.sra-instruments.com