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As seen on this graph, a start (via the corresponding GC contact or the START button on this page) launch the cycle at t_0 . The Modulation starts right ahead with a hot jet pulse in less than 2 ms. If you reset the SRA-Optimode with the START contact on, the cycle shall be running right now. Until t_1 , the first time coordinate of the mass-flow, is reached, the mass-flow is kept to the idle value. If $t_1 = 0$ s (strongly advised), the first mass-flow value is set at the same time as the first modulation. The mass-flow signal is linear with time between coordinates. When t_N , the last time coordinate, is reached, the system enter an "End" cycle, the mass-flow is kept to its last known value, and the modulation is stopped. A STOP(via the corresponding button or the GC contact on a shift down) is required to go back to the idle cycle.

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
Configuration

Hot Jet	MF Check Points	Idle Mass Flow	
Modulation (s) <input type="text" value="20"/>	<input type="text" value="2"/>	Mass Flow Idle (%) <input type="text" value="3"/>	<input type="button" value="Download"/>
HotJet (s) <input type="text" value="0.5"/>	<input type="text" value="00m00s500"/>		
Mass Flow			
Time (s) <input type="text" value="0"/>	<input type="text" value="00h15m00"/>		
<input type="text" value="70"/>	<input type="text" value="35"/>		
Mass Flow (%) <input type="text" value="70"/>	<input type="text" value="35"/>		

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Simple Run

From this page, you can set the parameters of a simple cycle and run it. The Optimode provide a 0V-24V pulse for hot jet control, and the aperture of an analogic mass-flow controller via a 4mA-20mA or 0-10V signal.

As seen on this graph, a start (via the corresponding GC contact or the START button on this page) launch the cycle at t_0 . The Modulation starts right ahead with a hot jet pulse in less than 2 ms. If you reset the SRA-Optimode with the START contact on, the cycle shall be running right now. Until t_1 , the first time coordinate of the mass-flow, is reached, the mass-flow is kept to the idle value. If $t_1 = 0$ s (strongly advised), the first mass-flow value is set at the same time as the first modulation. The mass-flow signal is linear with time between coordinates. When t_2 , the last time coordinate, is reached, the system enter an "End" cycle. the mass-flow is kept to its last known value, and the modulation is stopped. A STOP (via the corresponding button or the GC contact on a shift down) is required to go back to the idle Cycle.

The mass-flow and modulation parameters can be modified with the following formulary. Data won't be taken in account if the system is currently running. The number of mass-flow coordinate can be set between 1 and 10. The mass-flow time parameters must be sorted from smallest (0s advised) to longest. There is a 1 s time precision, up to approximately 5h. The mass-flow aperture is given in percentage, with a built in precision of 0.39%, and thus a 0.1% step in entry. The hot jet time, corresponding to the hot pulse, must be strictly superior to 0 ms, and the modulation time must be strictly superior to the hot jet. Both hot jet and modulation are defined with a 1 ms precision, up to a hundred hour.

If you enter an incorrect parameter, no value will be taken in account in the program. The current parameters appear above the entries.

- **Black** means the values are synchronized and your parameters have been accounted for.
- **Green** means the parameters haven't been updated due to a mistake somewhere in the formulary, and show you the current parameters used.
- **Red** also indicates parameters which haven't been taken in account, but are not part of the current cycle. (They are default values stored for further use)

All values are hard-saved only once the whole formulary is correct.

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Configuration

Hot Jet	MF Check Points	Idle Mass Flow	send
Modulation (s) <input type="text" value="20"/> 00m20s000	<input type="text" value="4"/>	Mass Flow Idle (%) <input type="text" value="3"/> 3	
HotJet (s) <input type="text" value="30"/> 00m00s500			
Mass Flow			
Time (s)	<input type="text" value="0"/> 00h00m00	<input type="text" value="900"/> 00h15m00	<input type="text" value="9060"/> 02h31m00
	<input type="text" value="50"/> 00h00m00	<input type="text" value="35"/> 00h00m00	<input type="text" value="5"/> 00h00m00
Mass Flow (%)	<input type="text" value="50"/>	<input type="text" value="35"/>	<input type="text" value="5"/>
	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text" value="0"/>

of any partial or total part of the software, hardware configuration, and htm material without SRA-Instruments agreemen

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Status : Stand By

Current Cycle
 Stand By
 Next Cycle
 Main Cycle
 GC Contact off
 Modulation off
 Mass Flow Value 3
 Elapsed Time
 Not Initialized

Start

Start ---

Overview

Next...

Welcome to the SRA-Optimode interface. The section in the left border allows navigation between the different sections.

- Overview (this section) offers a brief presentation of the configuration pages and of the status bar
- Simple-modulation enables simple runs of the Optimode with a fixed modulation period, configuration of the mass flow idle level, and mass-flow in-run control.
- Multi-modulation is used for more complex runs with different modulation and hot jet periods during a single run.
- Maintenance is a technical section where you can finely tune the internal Optimode dock, launch a defrost initiating cycle, test the mass-flow and hot-jet valve, and restore factory parameters
- IP-Config enables the setting of the board with fixed IP or automatically attributed DHCP address, and your network configuration.

Under the navigation tab is a status section, gathering the main information about the board and the software. The information displayed are quite straightforward, we shall nevertheless detail the "Current Cycle" and "Next Cycle" fields.

- Stand By is the default idle state of the device. The modulation is off and the mass-flow kept at idle value. The Idle Led is lit. Stand By is automatically entered at the reset of the board, if the "Start" connector is off (the connector status can be monitored on the "GC Start" status field. A down shift on the connector or the "STOP" button also engage Stand By mode.
- Main Cycle is a simple cycle with fixed modulation and mass-flow control. The Cycle Led is lit during operation. A Main Cycle shift to a Stand By mode if a stop (software button of connector) occurs, or to a "End" phase at the end of the cycle's duration. The time elapsed since the beginning of the run is displayed, along with the real time mass-flow value. The hot jet state (software) can also be monitored in real time.
- Special Cycle is similar to the Main Cycle, but enables different modulations in the same run. This cycle shall be launched from the proper web page (Multi-Modulation). If "Special Cycle" appears in the Next cycle field instead of "Main Cycle", reload the current web page with the "MAIN CYCLE" Button in the navigation bar.
- The End Cycle occurs at the end of an active cycle (Main or Multi). In this state, the modulation is off, and the mass-flow is kept at the last known value. To exit this state, one must call a stop (software or GC Start down shift), and go to Stand By mode. When the system enters the End Cycle, a "Press Stop" Message is displayed in the current Cycle field.
- Pause is a simple state entered with the corresponding software button, putting the device idle and on again. The Idle Led is lit during pause.
- Cold Up is a special function used at the beginning of your experiment. Via a brief sequence of hot pulses at 50% mass-flow aperture, the system is awakened and put in optimal conditions. Every hypothetical ice is removed by the hot jet, and the mass-flow arrival is cooled down. The cold up automatically ends in Stand By mode. The "Special Cycle" Led is lit when called, and "Initializing" is displayed.
- Config-MF is a special mode (the Special Led is lit) enabling fine tuning and test of the mass-flow controller and of the hot-jet. The Special Cycle Led is lit, this function is used at factory to set the 4 mA - 20 mA output trimmers.
- Factory Cycle is sometime selected. "Setting Factory Reset" is displayed, and the device parameters are set to their initial value. The IP configuration is not changed in between.

If something goes wrong in the software, an error message can occur. It shall be very unlikely, and you are advised to tell SRA Instrument of such occurrence. Using a software Start/Stop button a few times or resetting the board shall be enough to start anew.

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Status : Stand By

Current Cycle
 Stand By
 Next Cycle
 Main Cycle
 GC Contact off
 Modulation off
 Mass Flow Value 3
 Elapsed Time
 Not Initialized

Start

Simple Run

From this page, you can set the parameters of a simple cycle and run it. The Optimode provide a 0V-24V pulse for hot jet control, and the aperture of an analogic mass-flow controller via a 4mA-20mA or 0-10V signal.

As seen on this graph, a start (via the corresponding GC contact or the START button on this page) launch the cycle at t_0 . The Modulation starts right ahead with a hot jet pulse in less than 2 ms. If you reset the SRA-Optimode with the START contact on, the cycle shall be running right now. Until t_1 , the first time coordinate of the mass-flow, is reached, the mass-flow is kept to the idle value. If $t_1 = 0$ s (strongly advised), the first mass-flow value is set at the same time as the first modulation. The mass-flow signal is linear with time between coordinates. When t_n , the last time coordinate, is reached, the system enter an "End" cycle. the mass-flow is kept to its last known value, and the modulation is stopped. A STOP(via the corresponding button or the GC contact on a shift down) is required to go back to the idle Cycle.

The mass-flow and modulation parameters can be modified with the following formulary. Data won't be taken in account if the system is currently running. The number of mass-flow coordinate can be set between 1 and 10. The mass-flow time parameters must be sorted from smallest (0s advised) to longest. There is a 1 s time precision, up to approximately 5h. The mass-flow aperture is given in percentage, with a built in precision of 0.39%, and thus a 0.1% step in entry. The hot jet time, corresponding to the hot pulse, must be strictly superior to 0 ms, and the modulation time must be strictly superior to the hot jet. Both hot jet and modulation are defined with a 1 ms precision, up to a hundred hour.


if you enter an incorrect parameter, no value will be taken in account in the program. The current parameters appear above the entries.

- Black means the values are synchronized and your parameters have been accounted for.
- Green means the parameters haven't been updated due to a mistake somewhere in the formulary, and show you the current parameters used.
- Red also indicates parameters which haven't been taken in account, but are not part of the current cycle. (They are default values stored for further use)

All values are hard-saved only once the whole formulary is correct.

Configuration

Hot Jet	MF Check Points	Idle Mass Flow	
Modulation (s) <input style="width: 80px;" type="text" value="20"/>	<input style="width: 80px;" type="text" value="00m20s000"/> <input style="width: 40px;" type="text" value="2"/>	Mass Flow Idle (%) <input style="width: 80px;" type="text" value="3"/>	<input type="button" value="Download"/>
HotJet (s) <input style="width: 80px;" type="text" value="0.5"/>	<input style="width: 80px;" type="text" value="00m00s500"/>		
Mass Flow	<input style="width: 80px;" type="text" value="00h00m00"/> <input style="width: 80px;" type="text" value="00h15m00"/>		
Time (s) <input style="width: 80px;" type="text" value="0"/> <input style="width: 80px;" type="text" value="900"/>	<input style="width: 80px;" type="text" value="70"/> <input style="width: 80px;" type="text" value="35"/>		
Mass Flow (%) <input style="width: 80px;" type="text" value="70"/> <input style="width: 80px;" type="text" value="35"/>			


SRA Instruments 150, rue des Sources | 69280 Marcy l'Etoile | France
 Tel. +33 04 7844 2947 | Fax +33 04 7844 2962 | info@sra-instruments.com | www.sra-instruments.com
 Siege Social: 150, rue des Sources , 69280 Marcy l'Etoile
 SA a Directoire et a Conseil de Surveillance au capital de 150 000 Euros

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SRA Instruments
CHROMATOGRAPHIC SOLUTIONS

SRA Optimode

Release Version: 3.00
Compiling Date: Oct 1 2018 10:15:27

TCP/IP Version: v5.36
Serial Number: 2013

Menu

OVERVIEW

SIMPLE-MODULATION

MULTI-MODULATION

MAINTENANCE

IP CONFIGURATION

SPECIFICATIONS

Status : Stand By

Current Cycle
Stand By

Next Cycle

Multi Modulation

GC Contact off

Modulation off

Mass Flow Value 3

Elapsed Time

Not Initialized

Start

Start ---

Multi-Modulation

Multi-Modulation is a new feature in this Optimode. The modulation and Hot Jet are dynamically adjustable to your analysis. Like the mass-flow controller, you can enable up to 10 successive steps with different modulation parameters. The new modulation parameters are changed at the time associated. The modulation doesn't start directly anyway : the parameters are stored in memory, and become effective at the next automatic update which occurs at the end of each complete modulation. As a result, the different sections of a multi-modulation all share the same time shift. If you set the first modulation at time = 0 s, this offset is null, within the 2 ms precision.

The multi-modulation mode is similar to the simple modulation : it starts and stops in the same conditions, and the behavior of the mass-flow controller is the same. If the last mass-flow time parameter is smaller than the last modulation one, the mass-flow is kept at the same value until the cycle finish and enters "Ending" mode. If the mass-flow time is longer, the last modulation and pulse parameters are kept until the "Ending" Cycle. Also, the parameters are only updated and saved on hard memory once the whole formulary is correct.

The modulation can be disabled until next modulation time with the checkbox next to the time entry. When disabled, the modulation and pulse parameters are irrelevant, and are therefore not displayed. The next modulation shall start at the precise time set. When the Modulation is disabled in the program, no values (***) is displayed. However, the parameters are still available in memory for further use.

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

SRA Instruments
CHROMATOGRAPHIC SOLUTIONS

SRA Optimode

Release Version: 3.00 TCP/IP Version: v5.36
 Compiling Date: Oct 1 2018 10:15:27 Serial Number: 2013

Menu

OVERVIEW

SIMPLE-MODULATION

MULTI-MODULATION

MAINTENANCE

IP CONFIGURATION

SPECIFICATIONS

Status : Stand By

Current Cycle
Stand By
Next Cycle

Config MF
GC Contact off
Modulation off
Mass Flow Value 3
Elapsed Time
Not Initialized


Start

Start ---

Maintenance Section

The maintenance section is composed of 4 different functions.

- Mass Flow and Hot Jet Set is used in factory for testing and hardware configuration. The system enters a special mode when the formulary is submitted. The state of the hot jet and the mass-flow value can be configured at will. The mass-flow value is displayed in percentage, in 8-bits binary (the numerical value of the mass-flow), on a 0 - 10 V scale, and on a 4 - 20 mA scale, for direct control. The mass-flow value has no incidence on the idle value.
- Cold Up is a special cycle used to put the system in optimal conditions. The mass-flow is set to 50% aperture, and successive 1 s long modulation and hot jet pulse are triggered. We can therefore de-ice the cold jet, and cool the air entry. The parameters of this special cycle are factory written in hard memory, and can't be change.
- Timer shift enables the tuning of the SRA-Optimode internal dock. 1 ms is initially defined as 40.000 cycle on our timer, in accordance with a 80 MHz built-in quartz. As seen in the following figure, the span difference between the Optimode clock and your GC can lead to a difference in the modulation period, and to a shift in pulse over time. We can observe, by example, a 5 s shift over the course of 1h30. After characterizing the time shift and the total time, enter the required correction in the time shift formulary to synchronize your docks. The same unit (to your choice) is required, with relativ's integers only.
- The final module is used for memory save and recovery. Parameters are automatically saved after a correct entry, and this section enable factory recovering of different parameters : the timer definition of 1 ms, The idle level of the mass-flow, the Hot Jet parameters, and the mass-flow parameters.

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
SPECIFICATIONS

Status : Stand By

Current Cycle
Stand By
Next Cycle
Config MF
GC Contact off
Modulation off
Mass Flow Value 3
Elapsed Time
Not Initialized

IP Configuration

This page allows the configuration of the board's network settings. Be very careful, once you've submitted the new settings, the board may loose connectivity with the network. As soon as new settings are input, the board shall reset itself, resulting in the 3 status Led Blinking. Recovery options will be provided on the next page. The DHCP option is useful if you want the board to synchronize with your DHCP server. The AutoIP option (overrun by the DHCP) attributes a random IP to your board if it can't achieve connection.

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
SRA Optimode
Release Version: 3.00
 Compiling Date: Oct 1 2018 10:15:27

TCP/IP Version: v5.36
 Serial Number: 2013

Specifications

Hot Jet Signal (2D)

Specification Name	Specification Value
Hot Jet Modulation	
Modulation resolution (Mean)	0.1 µs
Hot Jet resolution (Mean)	0.1 µs
Modulation repeatability (Variance)	12.5 ns
Hot Jet repeatability (Variance)	12.5 ns
Modulation Step	1 ms
Hot Jet Step	1 ms
Modulation Length	2 ms to 1 h 40
Multi-Modulation steps	10
Multi-Modulation time parameter	up to 100 h
Initialisation delay at "Start" command	between 1 ms and 2 ms.
Mass Flow Controller	
Mass Flow Resolution	0.39 %
Mass Flow time parameter	up to 5h30 between each step up to 100h
Mass-Flow steps	10
Internal Clock	
maximum dock drift (after correction)	9 ms per hour.

Menu

- OVERVIEW
- SIMPLE-MODULATION
- MULTI-MODULATION
- MAINTENANCE
- IP CONFIGURATION
- SPECIFICATIONS

Status : Stand By

Current Cycle
 Stand By
 Next Cycle

Config MF
 GC Contact off
 Modulation off
 Mass Flow Value 3
 Elapsed Time
 Not Initialized

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