

BioVOC-2™

BioVOC-2 is a single-use device for cost-effectively collecting VOCs in breath onto sorbent tubes, for subsequent analysis by thermal desorption (TD).




Simple biological monitoring

Requiring minimal training, BioVOC-2 is ideal for any application where the presence of medically qualified staff would not be practical – such as large-scale studies of industrial workers or the general public.

Outstanding performance

BioVOC-2 has the same design as the original Bio-VOC™ (developed by the UK Health & Safety Executive), and uses a non-emitting, non-binding polymer in the sampler body.

This provides the same high performance as the previous model, minimising any risk of contaminant interference or adsorption of analytes.



The ability of BioVOC-2 to discriminate between known compounds is demonstrated in this analysis of the breath of a healthy volunteer. Performance is comparable to the previous Bio-VOC model.

Product highlights

- Collects the the last 129 mL of exhaled air – the most reliable indicator of VOC concentrations in the blood.
- Construction from optimised polymers minimises risk of analytical artefacts.
- Can be used multiple times for one participant within a single sampling session.
- Simple manual transfer of collected breath to TD sorbent tubes.



To collect a breath sample using BioVOC-2, the subject simply breathes into the disposable cardboard mouthpiece.

Applications of BioVOC-2

Breath sampling using BioVOC-2 is relevant to three main fields of study:

- **Clinical research:** Various medical conditions have been found to influence the profile of VOCs in the breath. Breath sampling using tools such as BioVOC-2 is therefore increasingly being recognised as suitable for large-scale screening of susceptible populations.
- **Occupational health monitoring:** Exposure to solvents and other workplace chemicals is easy to assess using BioVOC-2, and can be applied to both short-term and long-term monitoring.
- **Assessing exposure of the public to VOCs:** Chemicals in consumer products and construction materials can lead to significant and potentially harmful personal exposures, which can be monitored using BioVOC-2.

For more information on the use of BioVOC-2 and further details of typical applications, see [Application Note 147](#).

Click the link or scan the code to watch our [short video](#):



Ordering information

Part number	Description
C-BIO02-01	BioVOC-2 breath sampler*, pk 1
C-BIO02-10	BioVOC-2 breath sampler*, pk 10
C-BIO02-100	BioVOC-2 breath sampler*, pk 100
C-B100-M	Disposable cardboard mouthpieces, for BioVOC-2, pk 100
C2-CAXX-5149	Inert-coated sorbent tubes, 'Biomonitoring', pk 10

* Includes the corresponding number of cardboard mouthpieces.


Markes International

UK: Gwaun Elai Medi-Science Campus, Llantrisant, RCT, CF72 8XL
USA: 2355 Gold Meadow Way, Gold River, Sacramento, California 95670
Germany: Bieberer Straße 1-7, 63065 Offenbach am Main
P.R. China: No. 1 Building, No. 7 Guiqing Road, Xuhui District, Shanghai 200233
E: enquiries@markes.com **W:** www.markes.com

“Breath sampling using Markes’ Bio-VOC has been the cornerstone of our medical research for almost five years, thanks to its rapid, simple, and (most important) patient-friendly operation”

Professor James Covington

School of Engineering, University of Warwick



Markes International supplies a wide range of industry-standard 3½" x ¼" sorbent-packed tubes compatible with BioVOC-2.



Order your
BioVOC-2 today

Please contact your regional distributor or one of our offices (details below).

T: +44 (0)1443 230935
T: +1 866-483-5684 (toll-free)
T: +49 (0)69 6681089-10
T: +86 21 5465 1216

